



**MOBILITY DEVICES FOR PETS**

# PRODUCT GUIDE

DERRICK CAMPANA, CO (703) 474-6204



## ▣ PRODUCT GUIDE

- STIFLE BRACES
- CARPAL/TARSAL BRACES
- HOCK BRACES
- HIP BRACES
- PROSTHETICS
- CONTRACTURE BRACING
- DEFORMITY BRACING
- SOFT GOODS

## ▣ ORDERING INFO

- CALL: DERRICK CAMPANA, CO  
@ (703) 474-6204
- FAX ORDER FORM TO: (866) 885-8111
- E-Mail: [animalorthocare@hotmail.com](mailto:animalorthocare@hotmail.com)
- Address: 4200A Lafayette Center Dr.  
Chantilly, VA 20151
- Website: [www.animalorthocare.com](http://www.animalorthocare.com)
- PRICE GUIDE ON NEXT PAGE



# PRICE GUIDE

ALL PRICES INCLUDE  
CASTING KIT

CASTING KIT only- \$60 (kit is **FREE** with brace order and includes casting tape).

STIFLE-\$725 (any size or style)

CARPAL/TARSAL-\$545-\$595

HOCK- \$545-\$595

HIP- \$895

PROSTHETICS- \$1,250

CONTRACTURE- \$650-\$950

DEFORMITY-\$545-\$795

SOFT GOODS- PLEASE CALL

\* ALL joint types are equal in price minus ultraflex ® joints.

\* ALL PRICES ARE SUBJECT TO CHANGE AND DO NOT INCLUDE SHIPPING. Please call (703) 474-6204 for up-to-date pricing.

## ▣ Casting Kit Materials



# STIFLE BRACES

Used for CCL ruptures, post-op care, patellar issues and other instabilities.

STIFLE BRACE WITH "ROM" JOINT.

THIGH COMPONENT.

ADJUSTABLE "ROM" KNEE JOINT WITH SET SCREWS.

REMOVABLE SECONDARY SUSPENSION SYSTEM TO PREVENT MIGRATION.

STIFLE BRACE FOR MCL/LCL AND/OR PCL DEFICIENCIES.

TIBIAL AND POPLITEAL STRAPS TO PREVENT CRANIAL/CAUDAL DRAWER MOTIONS.

ALUMINUM UPRIGHTS TO PREVENT MEDIAL/LATERAL AND ROTATIONAL ISSUES.

STIFLE BRACE USED FOR PATELLAR LUXATION.

"ROM" JOINT TO LIMIT FLEXION/EXTENSION AND PREVENT PATELLAR DISLOCATION.



# CARPAL BRACES

Used for carpal hyperextension, varus/valgus deformities

"Carpal brace" with hinge and adjustable flexion strap. Strap can be elongated to allow for more motion/function, while preventing collapse of the carpus.

"check strap" to limit flexion.

"Tamarack" hinge to reinforce normal ligamentous structure.

CARPAL/PAW BRACE W/ROM

ADJUSTABLE RANGE OF MOTION JOINT.


"ROCKER BOTTOM" SOLE TO NORMALIZE GAIT PATTERN.

Carpal/Paw brace. Paw section added for maximum support. Can be removed to allow for increased function.

3-ply sock for added comfort/protection.

Adjustable check strap to control flexion and extension of hock.


"Tamarack" hinge.




# HOCK BRACES

Used for Achilles tendon ruptures or to replace external fixators.


"Hock Brace" with adjustable range of motion hinge to allow for set or limited motion. Hinge can be "locked out" or adjusted as desired with an included allen wrench. Mainly used for achilles tendon injuries or in replacement of an external fixator.



HEAVY DUTY CUSTOM HOCK/PAW BRACE WITH "ROM" JOINTS.



HOCK/PAW BRACE WITH "TAMARAC" JOINTS & FLEXION ADJUSTMENT TABS.



# HIP BRACES

Used for hip dysplasia and post operative support.



# PROSTHETICS


2-3 inches of radius/ulna for forelimb and as close to hock as possible for rear limb must remain in order to fabricate a custom prosthesis.





# CONTRACTURE BRACING

Ultraflex<sup>®</sup> joints are used to flex or extend contracted tendons.



# DEFORMITY BRACING



# SOFT GOODS

We can order and fit any product found online.

